

Product	Price	Ease of Use (Time) Harvest + Processing	Type of Graft	Viability of Graft	Addresses Sizing	Washing / Filtering	Centrifuge	All-Inclusive Kit	510K
intelliFat <sup>TM</sup> (Cellmyx)	\$895	30-40mins	Millifat (1.0mm)  SVT 200-500µm	Millifat: 90+% (No loss of final product)  SVT: 60-80% (less than 30% loss)	Yes	Yes for Millifat and SVT without centrifuging (like small Puregraft)	No	Yes	Yes
LipoGems (Stryker)	\$1,600	60-90mins	Microfat (400-600µm)	90+% (sacrifice 30cc to get 10cc final product)	Yes	Yes via saline wash	No	Yes	Yes
LipoPro (Ranfac)	\$550- \$750	45-60mins	Millifat	Unknown (sacrifice 30cc to get 10cc final product)	No	Yes via centrifugation	Yes	Yes	Yes
Emcyte	\$400	45-60mins	Millifat	Unknown	No	Yes via centrifugation	Yes	No	Yes
Harvest	\$500	45-60mins	Millifat	Unknown	No	Yes via centrifugation	Yes	Yes	Yes
Tulip	\$950	45-60mins	Millifat and Nanofat (dead fat)	Millifat: 60-80% (No loss of final product)  Nanofat: <20%	Yes	No unless using FatPress = extra \$\$\$	No	No	No